

# Bed-o'Cobs® 1/8" THE BEDDING OF CHOICE

## Keep Your Animals' Environment Dry With 1/8" Bed-o'Cobs

### GREATER SURFACE AREA

- Four times more surface area than 1/4" Bed-o'Cobs
- Will not cling to animals or cage
- Faster wicking absorption

### NATURAL ABSORPTION & ODOR CONTROL

Bed-o'Cobs is proven to encapsulate the urea salt crystals found in animal urine, controlling ammonia levels. Its sponge-like characteristics then encapsulate ammonia and trap associated odors.

### PERFECT FOR ISOLATOR AND IVC CAGES

Porous cob material allows liquids to flow to the bottom of the cage then absorbs liquids from the bottom up.

### ENVIRONMENTALLY FRIENDLY

This product is renewable, biodegradable, incinerable, compostable and is a product of annual crop production.

### OTHER SPECIFICATIONS

- No variables added since 1962
- Consistently the highest quality corn cob bedding
- Quarterly independent analysis testing

### PACKAGING:

- 1.25 cu ft / 35 L Bags
- Enough to fill 111 standard (78 sq in) shoebox cages 1/4" deep



P.O. Box 119 Maumee, Ohio 43537  
1-866-234-0505, Ext. 6511  
[www.AndersonsLabBedding.com](http://www.AndersonsLabBedding.com)

# FULL LINE OF BEDDING PRODUCTS FOR LABORATORY ANIMALS

|                        | DESCRIPTION                                              | PRODUCT CODE | CAGE FILLS | PACKAGING                          | PALLET STACK |
|------------------------|----------------------------------------------------------|--------------|------------|------------------------------------|--------------|
| BEDDING                | <b>Bed-o'Cobs® 1/4 in.</b> — 1/4" Cob*                   | 4B           | 111        | 1.25 cu ft/35 L bag                | 70           |
|                        | <b>Bed-o'Cobs® 1/8 in.</b> — 1/8" Cob                    | 8B           | 111        | 1.25 cu ft/35 L bag                | 60           |
|                        | <b>Bed-o'Cobs® Combo</b> — 1/8" & 1/4" Cob               | CB           | 111        | 1.25 cu ft/35 L bag                | 70           |
|                        | <b>Pure-o'Cel™</b> — Paper Squares*                      | PB           | 111        | 1.25 cu ft/35 L bag                | 80           |
|                        | <b>XTRAFine™</b> — Finely Ground Extra Absorbent Cob     | XFB          | 111        | 1.25 cu ft/35 L bag                | 60           |
| ENRICHMENT             | <b>Enrich-n'Nest™</b> — Paper Rolls                      | EB           | 1600       | 1.25 cu ft/35 L bag                | 80           |
|                        | <b>Crink-I'Nest™</b> — Crinkled Paper (kraft or white)   | CNK CNW      | 1600-2000  | 22 lb case                         | 15           |
| BEDDING AND ENRICHMENT | <b>Bed-o'Cobs/Pure-o'Cel™</b> — 1/4" Cob & Paper Squares | 4PB          | 111        | 1.25 cu ft/35 L bag                | 80           |
|                        | <b>Enrich-o'Cobs™</b> — 1/4" Cob & Paper Rolls           | 4EB          | 111        | 1.25 cu ft/35 L bag                | 80           |
|                        | <b>Enrich-n'Pure™</b> — Paper Rolls & Squares            | EPB          | 111        | 1.25 cu ft/35 L bag                | 80           |
| SPECIALTY              | <b>Surgical Suite Pad™</b> — Surgical Pad                | RPS1         | -          | 400/case                           | 25           |
|                        | <b>Recovery Pad™</b> — Post-Operative Pads               | RPM1<br>RPR1 | -          | Mouse (400/case)<br>Rat (250/case) | 25           |
|                        | <b>LabFresh™</b> — Extra Absorbent Litter                | LF20         | -          | 20 lb bag                          | 77           |
|                        | <b>Pel-o'Cobs®</b> — Pelleted Cob Bedding                | PELB         | -          | 1.25 cu ft/35 L bag                | 50           |
|                        | <b>Alf-a'Cobs®</b> — Pelleted Deodorizing Cob Bedding    | ALFB         | -          | 1.25 cu ft/35 L bag                | 50           |

\*Also available qualified

FOR THE NAME OF YOUR NEAREST DISTRIBUTOR - CALL US TODAY AT 1-866-234-0505 Ext. 6511



Check out our new website!

P.O. Box 119 Maumee, Ohio 43537  
1-866-234-0505, Ext. 6511  
[www.AndersonsLabBedding.com](http://www.AndersonsLabBedding.com)

