



Bed-o'Cobs[®] COMBO THE BEDDING OF CHOICE

Our Combination Bed-o'Cobs Delivers The Ideal Mix Of Comfort and Cleanliness

A BLEND OF OUR 1/8" AND 1/4" COBS THAT IS IDEAL FOR MICE, RATS, GUINEA PIGS AND OTHER SMALL ANIMALS

The blend of both particle sizes enhances enrichment while filling in the spaces and offering more surface area for fast absorption and ammonia control.

NATURAL ABSORPTION & ODOR CONTROL

Bed-o'Cobs is proven to encapsulate the urea salt crystals found in animal urine, controlling ammonia levels. Its sponge-like characteristics encapsulate ammonia and trap associated odors.

PERFECT FOR ISOLATOR AND IVC CAGES

Porous cob material allows liquids to flow to the bottom of the cage then absorbs liquids from the bottom up.

ENVIRONMENTALLY FRIENDLY

This product is renewable, biodegradable, incinerable, compostable and is a product of annual crop production.

OTHER SPECIFICATIONS

- No variables added since 1962
- Consistently the highest quality corn cob bedding
- Quarterly independent analysis testing

PACKAGING:

- 1.25 cu ft / 35 L Bags
- Enough to fill 111 standard (78 sq in) shoebox cages 1/4" deep



P.O. Box 119 Maumee, Ohio 43537
1-866-234-0505, Ext. 6511
www.AndersonsLabBedding.com

FULL LINE OF BEDDING PRODUCTS FOR LABORATORY ANIMALS

	DESCRIPTION	PRODUCT CODE	CAGE FILLS	PACKAGING	PALLET STACK
BEDDING	Bed-o'Cobs® 1/4 in. — 1/4" Cob*	4B	111	1.25 cu ft/35 L bag	70
	Bed-o'Cobs® 1/8 in. — 1/8" Cob	8B	111	1.25 cu ft/35 L bag	60
	Bed-o'Cobs® Combo — 1/8" & 1/4" Cob	CB	111	1.25 cu ft/35 L bag	70
	Pure-o'Cel™ — Paper Squares*	PB	111	1.25 cu ft/35 L bag	80
	XTRAFine™ — Finely Ground Extra Absorbent Cob	XFB	111	1.25 cu ft/35 L bag	60
ENRICHMENT	Enrich-n'Nest™ — Paper Rolls	EB	1600	1.25 cu ft/35 L bag	80
	Crink-I'Nest™ — Crinkled Paper (kraft or white)	CNK CNW	1600-2000	22 lb case	15
BEDDING AND ENRICHMENT	Bed-o'Cobs/Pure-o'Cel™ — 1/4" Cob & Paper Squares	4PB	111	1.25 cu ft/35 L bag	80
	Enrich-o'Cobs™ — 1/4" Cob & Paper Rolls	4EB	111	1.25 cu ft/35 L bag	80
	Enrich-n'Pure™ — Paper Rolls & Squares	EPB	111	1.25 cu ft/35 L bag	80
SPECIALTY	Surgical Suite Pad™ — Surgical Pad	RPS1	-	400/case	25
	Recovery Pad™ — Post-Operative Pads	RPM1 RPR1	-	Mouse (400/case) Rat (250/case)	25
	LabFresh™ — Extra Absorbent Litter	LF20	-	20 lb bag	77
	Pel-o'Cobs® — Pelleted Cob Bedding	PELB	-	1.25 cu ft/35 L bag	50
	Alf-a'Cobs® — Pelleted Deodorizing Cob Bedding	ALFB	-	1.25 cu ft/35 L bag	50

*Also available qualified

FOR THE NAME OF YOUR NEAREST DISTRIBUTOR - CALL US TODAY AT 1-866-234-0505 Ext. 6511



Check out our new website!

P.O. Box 119 Maumee, Ohio 43537
1-866-234-0505, Ext. 6511
www.AndersonsLabBedding.com

